



Stereo Microscope TMS-1420



\$650.00

Microscopes	Stereo Microscope TMS-1000 Series	Stereo Microscope TMS-1400 Series	Stereo Microscopes
-------------	-----------------------------------	-----------------------------------	--------------------

Stereo Microscope TMS-1420 suitable for school education, circuit board inspection, and repair, surface mount technology work, electronics inspection, coin collecting, gemology and gemstone setting, engraving, inspection and repair of small parts. The binocular viewing head and V-shaped objective. The stand has the same ergonomic design with both the incident and transmitted lighting system are rheostat controlled by halogen lamp and Incandescent lamp that allows the illumination to be adjusted to suit the sample type. Provide upright, un-reversed 3D images with high resolution and high quality.

[Add To Cart](#)

Product Description

Stereo Microscope TMS-1420 suitable for school education, circuit board inspection, and repair, surface mount technology work, electronics inspection, coin collecting, gemology and gemstone setting, engraving, inspection and repair of small parts. The binocular viewing head and V-shaped objective. The stand has the same ergonomic design with both the incident and transmitted lighting system are rheostat controlled by halogen lamp and Incandescent lamp that allows the illumination to be adjusted to suit the sample type. Provide upright, un-reversed 3D images with high resolution and high quality.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now