



Tabletop Autoclave TA-314



\$3,500.00

Steam Autoclave TA-300 Series	Steam Sterilisers & Autoclaves	Tabletop Front-loading Autoclaves
-------------------------------	--------------------------------	-----------------------------------

Tabletop Autoclave TA-314 front loading 'Benchtop' autoclave range is available in three sizes. It require no plumbing, and only a standard single phase power supply. Build in LCD display with 8 programs automatically feeding water and draining water system which is easy to operate and complete automatic operation, check and monitor by computer.

Add To Cart

Product Description

Tabletop Autoclave TA-314 front loading 'Benchtop' autoclave range is available in three sizes. It require no plumbing, and only a standard single phase power supply. Build in LCD display with 8 programs automatically feeding water and draining water system which is easy to operate and complete automatic operation, check and monitor by computer.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now