



Ultra-Low Chest Freezers -152°C ULTF-4258



\$26,700.00

Laboratory Refrigerators and Freezers	Ultra-Low Chest Freezers -152°C ULTF-4000	Ultra-Low Chest Freezers -86°C/-152°C/-164°C	Ultra-Low Freezers	Ultra-Low Freezers -86°C/-152°C/-164°C
---------------------------------------	---	--	--------------------	--

Ultra-Low Chest Freezers -152°C ULTF-4000 with VIP vacuum three-dimensional thermal insulation system and one pole lubrication compressor: refrigeration technology. Secure sample at the ultra-low and cryogenic conditions at -152°C. Equipped with many alarms function Visual & Audio high and low-temperature alarm, High/low voltage, sensor failure alarm, Low battery, Condenser heat dissipation, System failure. Features with single-pole lubricated compressor refrigeration technology for better-quality refrigeration capacity. Allow maximum capacity up to 258 L with a temperature range of -110°C to -152°C by offering excellent temperature uniformity and stability. Available with Upright racks, alarm lamp, voltage compensation, remote alarm, and liquid N2 backup system as optional accessories

[Add To Cart](#)

Product Description

Ultra-Low Chest Freezers -152°C ULTF-4258 with VIP vacuum three-dimensional thermal insulation system and one pole lubrication compressor: refrigeration technology. Secure sample at the ultra-low and cryogenic conditions at -152°C. Equipped with many alarms function Visual & Audio high and low-temperature alarm, High/low voltage, sensor failure alarm, Low battery, Condenser heat dissipation, System failure. Features with single-pole lubricated compressor refrigeration technology for better-quality refrigeration capacity. Allow maximum capacity up to 258 L with a temperature range of -110°C to -152°C by offering excellent temperature uniformity and stability. Available with Upright racks, alarm lamp, voltage compensation, remote alarm, and liquid N2 backup system as optional accessories

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now