



Ultra-Low Chest Freezers -164°C ULTF-5128



\$27,300.00

Laboratory Refrigerators and Freezers	Ultra-Low Chest Freezers -164°C ULTF-5000	Ultra-Low Chest Freezers -86°C/-152°C/-164°C	Ultra-Low Freezers	Ultra-Low Freezers -86°C/-152°C/-164°C
---------------------------------------	---	--	--------------------	--

Ultra-Low Chest Freezers -164°C ULTF-5128 build-in foaming technology, and VIP insulation which greatly improves the insulation effect, VIP vacuum insulation board, tightly locks the cold air to ensure the great insulation effect. Features with single-pole lubricated compressor refrigeration technology for better-quality refrigeration capacity. Allow maximum capacity up to 128 L with a temperature range of -110°C to -164°C by offering excellent temperature uniformity and stability. Available with Upright racks, alarm lamp, voltage compensation, remote alarm, and liquid N2 backup system as optional accessories.

[Add To Cart](#)

Product Description

Ultra-Low Chest Freezers -164°C ULTF-5128 build-in foaming technology, and VIP insulation which greatly improves the insulation effect, VIP vacuum insulation board, tightly locks the cold air to ensure the great insulation effect. Features with single-pole lubricated compressor refrigeration technology for better-quality refrigeration capacity. Allow maximum capacity up to 128 L with a temperature range of -110°C to -164°C by offering excellent temperature uniformity and stability. Available with Upright racks, alarm lamp, voltage compensation, remote alarm, and liquid N2 backup system as optional accessories.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now