



Ultra-Low Upright Freezers -86°C ULTF-2650



\$11,388.00

Laboratory Refrigerators and Freezers	Ultra-Low Freezers	Ultra-Low Freezers -86°C/-152°C/-164°C	Ultra-Low Upright Freezers -86°C ULTF-2000 Series	Ultra-Low Upright Freezers -86°C/-152°C/-164°C
---------------------------------------	--------------------	--	---	--

Ultra-Low Upright Freezers -86°C ULTF-2650 design with a cascade refrigeration system that uses two kinds of refrigerants having different boiling points, which run through their own independent freezing cycle and are joined by a heat exchanger. Features with filter for easy cleaning and single hand openable, Gasket heating, leveling caster, etc. Available with maximum capacity up to 650 L with a temperature range of -50 to-86°C.

[Add To Cart](#)

Product Description

Ultra-Low Upright Freezers -86°C ULTF-2650 design with a cascade refrigeration system that uses two kinds of refrigerants having different boiling points, which run through their own independent freezing cycle and are joined by a heat exchanger. Features with filter for easy cleaning and single hand openable, Gasket heating, leveling caster, etc. Available with maximum capacity up to 650 L with a temperature range of -50 to-86°C.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now