



Ultra-Low Upright Freezers

-86°C ULTF-3000(110V/60Hz)



\$12,300.00

Laboratory Refrigerators and Freezers	Ultra-Low Chest Freezers -86°C ULTF-3000	Ultra-Low Chest Freezers -86°C/-152°C/-164°C	Ultra-Low Freezers	Ultra-Low Freezers -86°C/-152°C/-164°C
---------------------------------------	--	--	--------------------	--

Ultra-Low Upright Freezers -86°C ULTF-3000 suitable for blood stations, hospitals, health and epidemic prevention system, scientific research institutes, biological engineering, a laboratory in colleges and universities, military-industrial enterprises, etc. Available with a maximum capacity up to 668 L. Microcomputer temperature control system to set the inner temperature at the range of -40 to -86°C.

[Add To Cart](#)

Product Description

Ultra-Low Upright Freezers -86°C ULTF-3000 suitable for blood stations, hospitals, health and epidemic prevention system, scientific research institutes, biological engineering, a laboratory in colleges and universities, military-industrial enterprises, etc. Available with a maximum capacity up to 668 L. Microcomputer temperature control system to set the inner temperature at the range of -40 to -86°C.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now