



## Upright Biomedical Freezers -25°C MF-100(110V/60Hz)



\$2,250.00

Laboratory Refrigerators and Freezers	Ultra-Low Freezers	Upright Biomedical Freezers -25°C MF-100	Upright Freezers -25°C
---------------------------------------	--------------------	--	------------------------

Upright Biomedical Freezers -25°C MF-100 is a single-door Upright freezer with high energy and temperature precision technology available with a maximum capacity of up to 270L. Built-in doors sealing, dust-proof, easy to clean functions with maximum storage space. The system materials are ABS engineering plastic liners with adjustable front screws that facilitate easy movement. Ideal for research applications, storage of low-temperature materials, and cryopreservation of blood plasma, vaccines, and other biological products. These biomedical freezers provide stable temperatures from -10 to -25°C. Features include digital temperature display, temperature, microcomputer temperature control, compressor and condenser, 7 drawers, thermal insulating foamed door, alarm system, and password protections that provide the most reliable preservation system.

[Add To Cart](#)

### Product Description

Upright Biomedical Freezers -25°C MF-100 is a single-door Upright freezer with high energy and temperature precision technology available with a maximum capacity of up to 270L. Built-in doors sealing, dust-proof, easy to clean functions with maximum storage space. The system materials are ABS engineering plastic liners with adjustable front screws that facilitate easy movement. Ideal for research applications, storage of low-temperature materials, and cryopreservation of blood plasma, vaccines, and other biological products. These biomedical freezers provide stable temperatures from -10 to -25°C. Features include digital temperature display, temperature, microcomputer temperature control, compressor and condenser, 7 drawers, thermal insulating foamed door, alarm system, and password protections that provide the most reliable preservation system.

### Similar Products



Buy Now



Buy Now



Buy Now



Buy Now