



Upright Biomedical Freezers -25°C MF-260



\$3,201.00

Laboratory Refrigerators and Freezers	Medical Freezers	Ultra-Low Freezers	Upright Biomedical Freezers -25°C MF-200
---------------------------------------	------------------	--------------------	--

Upright Biomedical Freezers -25°C MF-260 have been designed to maintain a temperature range of -10°C to -25°C. Single door low-temperature freezer as they contain shelves and draws, letting foods to be quickly sorted. Available with maximum capacity up to 263L. Equipped with automatic magnetic sealed door closes to ensure and heightens the refrigeration performance and operational reliability, and password protection to control operating parameters from being modified, Robust and durable framing, antibacterial powdered inner and outer surface, refrigeration system uses highly efficient and eco-friendly and fluorine-free refrigerant. Rechargeable battery, which can display the main information such as alarm and the internal temperature of the cabinet when the main power is cut off and mechanical lock to ensure safety.

Add To Cart

Product Description

Upright Biomedical Freezers -25°C MF-260 have been designed to maintain a temperature range of -10°C to -25°C. Single door low-temperature freezer as they contain shelves and draws, letting foods to be quickly sorted. Available with maximum capacity up to 263L. Equipped with automatic magnetic sealed door closes to ensure and heightens the refrigeration performance and operational reliability, and password protection to control operating parameters from being modified, Robust and durable framing, antibacterial powdered inner and outer surface, refrigeration system uses highly efficient and eco-friendly and fluorine-free refrigerant. Rechargeable battery, which can display the main information such as alarm and the internal temperature of the cabinet when the main power is cut off and mechanical lock to ensure safety.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now