



Upright Deep Freezers -40°C DF-260



\$3,348.00

| | | | | |
|------------------------|---|-----------------------|--------------------------------|---|
| Deep Freezers -40°C | Laboratory Refrigerators and Freezers | Ultra-Low Freezers | Upright Deep Freezers -40°C | Upright Deep Freezers -40°C DF-200 Series |
|------------------------|---|-----------------------|--------------------------------|---|

Upright Deep Freezers -40°C DF-260 is well organized single-door low-temperature freezer as they contain shelves and draws, letting foods be quickly sorted. The maximum capacities up to 263L with a temperature range of -30°C to -40°C. Equipped with automatic magnetic sealed door closes to ensure and heightens the refrigeration performance and operational reliability, and password protection to control operating parameters from being modified, robust and durable framing, antibacterial powdered inner and outer surface, refrigeration system uses highly efficient and eco-friendly and fluorine-free refrigerant. Furnished with a rechargeable battery, which can display the main information such as alarm and the internal temperature of the cabinet when the main power is cut off Built-in mechanical lock to ensure safety.

Add To Cart

Product Description

Upright Deep Freezers -40°C DF-260 is well organized single-door low-temperature freezer as they contain shelves and draws, letting foods be quickly sorted. The maximum capacities up to 263L with a temperature range of -30°C to -40°C. Equipped with automatic magnetic sealed door closes to ensure and heightens the refrigeration performance and operational reliability, and password protection to control operating parameters from being modified, robust and durable framing, antibacterial powdered inner and outer surface, refrigeration system uses highly efficient and eco-friendly and fluorine-free refrigerant. Furnished with a rechargeable battery, which can display the main information such as alarm and the internal temperature of the cabinet when the main power is cut off Built-in mechanical lock to ensure safety.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now