

UV-VIS Spectro DS-1200



\$9,740.00



Double Beam UV-VIS Spectrophotometer	Double Beam UV-VIS Spectrophotometer DS-1000 Series	Spectrophotometer	UV-VIS Spectrophotometers
--------------------------------------	---	-------------------	---------------------------

Double Beam UV-VIS Spectrophotometer DS-1200 offer significantly better-quality and usability, a user-friendly touch screen display, design with 1200-line grating double beam spectrophotometer to ensure high resolution with wavelength availability of 190 to 1100 nm and 1.8 nm bandwidth, D2 Lamp & W Lamp. Useful function is photometric measurement, quantitative measurement, spectral scanning, dynamics, DNA/ protein testing, multi-wavelength testing and data printing functions.

[Add To Cart](#)

Product Description

Double Beam UV-VIS Spectrophotometer DS-1200 offer significantly better-quality and usability, a user-friendly touch screen display, design with 1200-line grating double beam spectrophotometer to ensure high resolution with wavelength availability of 190 to 1100 nm and 1.8 nm bandwidth, D2 Lamp & W Lamp. Useful function is photometric measurement, quantitative measurement, spectral scanning, dynamics, DNA/ protein testing, multi-wavelength testing and data printing functions.

Similar Products



[Buy Now](#)



[Buy Now](#)



[Buy Now](#)



[Buy Now](#)