

UV-VIS Spectro US-1600-220V



\$10,250.00



Single Beam UV-VIS Spectrophotometer US-1000 Series	Single Beam UV-VIS Spectrophotometers	Spectrophotometer	UV-VIS Spectrophotometers
---	---------------------------------------	-------------------	---------------------------

Single Beam UV-VIS Spectrophotometer US-1600 offers superb application functions with accuracy and precision. It is the best instrument for various applications. With 1200 line grating single beam to ensure high resolution with wavelength availability of 190 to 1000 nm and 4 nm bandwidth, tungsten, and D2 lamp. A useful function is an absorption, transmittance, and wavelength with automatic switching mode.

Add To Cart

Product Description

Single Beam UV-VIS Spectrophotometer US-1600 offers superb application functions with accuracy and precision. It is the best instrument for various applications. With 1200 line grating single beam to ensure high resolution with wavelength availability of 190 to 1000 nm and 4 nm bandwidth, tungsten, and D2 lamp. A useful function is an absorption, transmittance, and wavelength with automatic switching mode.


Similar Products



Buy Now



Buy Now



Buy Now



Buy Now