

UV-VIS Spectro US-2300-110 V



\$7,551.00



Single Beam UV-VIS Spectrophotometer US-2000 Series	Single Beam UV-VIS Spectrophotometers	Spectrophotometer	UV-VIS Spectrophotometers
-----------------------------------------------------	---------------------------------------	-------------------	---------------------------

Single Beam UV-VIS Spectrophotometer US-2300 provide sufficient light source for exact measurement. Our line of spectrophotometers is intend to fit into for QC, common research, pharmaceuticals, biochemical and clinical laboratory. Build-in halogen tungsten light source and silicon photodiode detector. We offer soft key pad, with continuous wavelength ranging from 190 to 1100 nm with variable spectral bandwidth 0.5, 1, 2, 4 and 5 nm and 160 to 200 mm focal length and used to measure light absorbance of chemical substance. This series is suitable for measuring photometry, quantity, multi-wavelength, spectrum, kinetics.

[Add To Cart](#)

Product Description

Single Beam UV-VIS Spectrophotometer US-2300 provide sufficient light source for exact measurement. Our line of spectrophotometers is intend to fit into for QC, common research, pharmaceuticals, biochemical and clinical laboratory. Build-in halogen tungsten light source and silicon photodiode detector. We offer soft key pad, with continuous wavelength ranging from 190 to 1100 nm with variable spectral bandwidth 0.5, 1, 2, 4 and 5 nm and 160 to 200 mm focal length and used to measure light absorbance of chemical substance. This series is suitable for measuring photometry, quantity, multi-wavelength, spectrum, kinetics.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now