



## Vertical Autoclave/Steam sterilizers VA-1800



\$7,970.00

Steam  
Sterilisers &  
Autoclaves

Vertical  
Autoclave

Vertical  
Autoclave/Steam  
sterilizers  
VA-1000 Series

Vertical Autoclave/Steam sterilizers VA-1800 floor standing self-contained sterilizer units that require no plumbing and operate from a standard single-phase power supply. Equipped with an effective air-cooling system for liquid sterilisation. This speeds cooling time by up to 40% compared to standard models. The Texas Vertical Autoclaves are available from 36 litres, 54 litres, 80 litres, and 100 litres chamber capacity. All are built with castor wheels as standard making them easy to move. Widely used for sterilisation applications in laboratory, hospitals, and clinics as well as for industrial sterilisation applications in biotechnology and pharmaceutical industries. Comes with sterilization cycle including solid cycle, Liquid cycle, Liquid with warming cycle, Agar cycle, Solid with drying (for D model), Customized cycle.

[Add To Cart](#)

### Product Description

Vertical Autoclave/Steam sterilizers VA-1800 floor standing self-contained sterilizer units that require no plumbing and operate from a standard single-phase power supply. Equipped with an effective air-cooling system for liquid sterilisation. This speeds cooling time by up to 40% compared to standard models. The Texas Vertical Autoclaves are available from 36 litres, 54 litres, 80 litres, and 100 litres chamber capacity. All are built with castor wheels as standard making them easy to move. Widely used for sterilisation applications in laboratory, hospitals, and clinics as well as for industrial sterilisation applications in biotechnology and pharmaceutical industries. Comes with sterilization cycle including solid cycle, Liquid cycle, Liquid with warming cycle, Agar cycle, Solid with drying (for D model), Customized cycle.

### Similar Products



Buy Now



Buy Now



Buy Now



Buy Now