



Vertical Autoclave/Steam sterilizers VA-2200



\$7,535.00

Steam
Sterilisers &
Autoclaves

Vertical
Autoclave

Vertical
Autoclave/Steam
sterilizers
VA-2000 Series

Vertical Autoclave/Steam sterilizers VA-2200 floor standing self-contained sterilizer units that require no plumbing and operate from a standard single-phase power supply. Equipped with an effective air-cooling system for liquid sterilization. This speeds cooling time by up to 40% compared to standard models. The Texas Vertical Autoclaves are available in 60 L chamber capacity. All are built with castor wheels as standard making them easy to move. Widely used for sterilization applications in laboratory, hospitals, and clinics as well as for industrial sterilization applications in biotechnology and pharmaceutical industries. Comes with sterilization cycle including Solid cycle, Liquid cycle, Liquid with warming cycle, Agar cycle, Solid with drying (for WC model), Medical waste cycle, Customized cycle.

[Add To Cart](#)

Product Description

Vertical Autoclave/Steam sterilizers VA-2200 floor standing self-contained sterilizer units that require no plumbing and operate from a standard single-phase power supply. Equipped with an effective air-cooling system for liquid sterilisation. This speeds cooling time by up to 40% compared to standard models. The Texas Vertical Autoclaves are available in 60 L chamber capacity. All are built with castor wheels as standard making them easy to move. Widely used for sterilization applications in laboratory, hospitals, and clinics as well as for industrial sterilisation applications in biotechnology and pharmaceutical industries. Comes with sterilization cycle including Solid cycle, Liquid cycle, Liquid with warming cycle, Agar cycle, Solid with drying (for WC model), Medical waste cycle, Customized cycle.

Similar Products

