



Vertical Autoclave/Steam sterilizers VA-2400



\$8,540.00

Steam Sterilisers & Autoclaves	Vertical Autoclave	Vertical Autoclave/Steam sterilizers VA-2000 Series
--------------------------------	--------------------	---

Vertical Autoclave/Steam sterilizers VA-2400 floor standing self-contained sterilizer units that require no plumbing and operate from a standard single-phase power supply. Equipped with an effective air-cooling system for liquid sterilization. This speeds cooling time by up to 40% compared to standard models. The Texas Vertical Autoclaves are available in 85 L chamber capacity. All are built with castor wheels as standard making them easy to move. Widely used for sterilization applications in laboratory, hospitals, and clinics as well as for industrial sterilization applications in biotechnology and pharmaceutical industries. Comes with sterilization cycle including Solid cycle, Liquid cycle, Liquid with warming cycle, Agar cycle, Solid with drying (for WC model), Medical waste cycle, Customized cycle.

Add To Cart

Product Description

Vertical Autoclave/Steam sterilizers VA-2400 floor standing self-contained sterilizer units that require no plumbing and operate from a standard single-phase power supply. Equipped with an effective air-cooling system for liquid sterilization. This speeds cooling time by up to 40% compared to standard models. The Texas Vertical Autoclaves are available in 85 L chamber capacity. All are built with castor wheels as standard making them easy to move. Widely used for sterilization applications in laboratory, hospitals, and clinics as well as for industrial sterilization applications in biotechnology and pharmaceutical industries. Comes with sterilization cycle including Solid cycle, Liquid cycle, Liquid with warming cycle, Agar cycle, Solid with drying (for WC model), Medical waste cycle, Customized cycle.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now