



Vertical Autoclave/Steam sterilizers VA-5100



\$7,841.00

Steam
Sterilisers &
Autoclaves

Vertical
Autoclave

Vertical
Autoclave/Steam
sterilizers
VA-5000 Series

Vertical Autoclave/Steam sterilizers VA-5100 floor standing self-contained sterilizer units that require no plumbing and operate from a standard single-phase power supply. Equipped with an effective air-cooling system for liquid sterilization. This speeds cooling time by up to 40% compared to standard models. The Texas Vertical Autoclaves is available from 100 L chamber capacity. All are built with castor wheels as standard making them easy to move. Widely used for sterilization applications in laboratory, hospitals, and clinics as well as for industrial sterilization applications in biotechnology and pharmaceutical industries.

[Add To Cart](#)

Product Description

Vertical Autoclave/Steam sterilizers VA-5100 floor standing self-contained sterilizer units that require no plumbing and operate from a standard single-phase power supply. Equipped with an effective air-cooling system for liquid sterilization. This speeds cooling time by up to 40% compared to standard models. The Texas Vertical Autoclaves is available from 100 L chamber capacity. All are built with castor wheels as standard making them easy to move. Widely used for sterilization applications in laboratory, hospitals, and clinics as well as for industrial sterilization applications in biotechnology and pharmaceutical industries.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now