



Vertical Laminar Flow Cabinet VLF-150

★★★★★

\$5,588.00

Laminar Flow Cabinets	Vertical Laminar Flow Cabinet VLF-100 Series	Vertical Laminar Flow Cabinets	Vertical Laminar Flow Cabinets VLF Series
-----------------------	--	--------------------------------	---

Model

Type

Cleanliness

Velocity

Material

Horizontal

Laminar

Flow

Cabinets

Series HLF

Series

Noise

Illumination

Power Consumption

Working Area

Overall Dimension

Weight

Electrical Requirement



Vertical Laminar Flow Cabinet VLF-150 provide hygienic and disinfected conditions for the processed material, recommended for a non-hazard sample. Available in three-single vertical operation mode used in medicine, pharmacy, and biological, chemical experiment, operating room, foodstuff, drinks, and cosmetic, stock rising, and aviation work.

Add To Cart

Product Description

Vertical Laminar Flow Cabinet VLF-150 provide hygienic and disinfected conditions for the processed material, recommended for a non-hazard sample. Available in three-single vertical operation mode used in medicine, pharmacy, and biological, chemical experiment, operating room, foodstuff, drinks, and cosmetic, stock rising, and aviation work. It is a general-purpose workbench, that delivers clean airflow from the top to the bottom in the vertical direction. Accommodates maximum no. of colony rate: ≤ 0.5 per utensil hour (90 mm utensil) and different air velocity ranging from 0.25-0.45m/s to 0.36-0.54m/s. Features with HEPA Filter, UV lamp, sliding glass door, ventilating slot, color plate, and LED source. Deliver best and economic free environment for inflexible contamination requirements on pharmaceutical, food industries, biotech, and much more Industries. Our flow cabinet enclosed on the sides and kept under constant pressure to prevent contamination.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now