



Vertical Laminar Flow Cabinets VLF-110

★★★★★

\$2,120.00

Laminar Flow Cabinets	Vertical Laminar Flow Cabinet VLF-100 Series	Vertical Laminar Flow Cabinets	Vertical Laminar Flow Cabinets VLF Series
-----------------------	--	--------------------------------	---

Model

Type

Cleanliness

Velocity

Material

Horizontal

Laminar

Flow

Cabinet

s HLF

Series

No. of colony

Noise

Illumination

Power Consumption

Vertical

Laminar

Flow

Cabinet

s
VLK-10
0
Series
Working Area
Micropl
ate
Sealers
Overall Dimension
Weight
Electrical Requirement
Catalogue No.

Vertical Laminar Flow Cabinet VLF-110 provide hygienic and disinfected conditions for the processed material, recommended for a non-hazard sample. Available in single-single vertical operation mode used in medicine, pharmacy, and biological, chemical experiment, operating room, foodstuff, drinks, and cosmetic, stock rising, and aviation work.

[Add To Cart](#)

Product Description

Vertical Laminar Flow Cabinet VLF-110 provide hygienic and disinfected conditions for the processed material, recommended for a non-hazard sample. Available in single-single vertical operation mode used in medicine, pharmacy, and biological, chemical experiment, operating room, foodstuff, drinks, and cosmetic, stock rising, and aviation work. It is a general-purpose workbench, that delivers clean airflow from the top to the bottom in the vertical direction. Accommodates maximum no. of colony rate: ≤ 0.5 per utensil hour (90 mm utensil) and different air velocity ranging from 0.25-0.45m/s to 0.36-0.54m/s. Features with HEPA Filter, UV lamp, sliding glass door, ventilating slot, color plate, and LED source. Deliver best and economic free environment for inflexible contamination requirements on pharmaceutical, food industries, biotech, and much more Industries. Our flow cabinet enclosed on the sides and kept under constant pressure to prevent contamination.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now