



Wide Field Zoom Stereo Microscope TMS-3265



\$3,800.00

Advanced Stereo Zoom Microscope TMS-3000 Series	Microscopes	Stereo Microscopes	Wide Field Zoom Stereo Microscope TMS-3200 Series
---	-------------	--------------------	---

Wide Field Zoom Stereo Microscope TMS-3265 designed for large viewing field without eyepiece for sharp and stereoscopic images from the lens in microscope head and has excellent objective to choose 2X, 4X, 6X, with optional choice 8X, 10X, 15X and 6XSL objectives. TMS-3200 Series is capable of providing strong images with a broad field of view, the multipurpose and perfect for basic studies of sample properties. Come with LED lights, allowing sample to be illuminated for a more comprehensive field of view. Suitable for Industrial, Medical and Scientific Research Areas, circuit board repair, circuit board inspection, surface mount technology work, electronics inspection, coin collecting, gemology and gemstone setting, engraving, and repair and inspection of small parts. Provide upright, un-reversed 3 D images. Used to examine the huge specimen under 4X-12X zoom magnification objective.

[Add To Cart](#)

Product Description

Wide Field Zoom Stereo Microscope TMS-3265 designed for large viewing field without eyepiece for sharp and stereoscopic images from the lens in microscope head and has excellent objective to choose 2X, 4X, 6X, with optional choice 8X, 10X, 15X and 6XSL objectives. TMS-3200 Series is capable of providing strong images with a broad field of view, the multipurpose and perfect for basic studies of sample properties. Come with LED lights, allowing sample to be illuminated for a more comprehensive field of view. Suitable for Industrial, Medical and Scientific Research Areas, circuit board repair, circuit board inspection, surface mount technology work, electronics inspection, coin collecting, gemology and gemstone setting, engraving, and repair and inspection of small parts. Provide upright, un-reversed 3 D images. Used to examine the huge specimen under 4X-12X zoom magnification objective.

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now