



WP-50(220V/50Hz)



\$1,680.00

Baths &
Circulators

Circulating
Water Bath
WP Series

Water Baths

1. LED Temperature display, resolution 0.1°C
2. Temperature control of the external system
3. PID temperature control
4. Low noise during the operation
5. Convenient to use on a table
6. Cooling coil (Option)
7. Adjustable high-temperature cut-off by heater and bath temperature sensor
8. Dry running protection
9. Protection from excess electric current
10. Water level sensor (Option)

[Add To Cart](#)

Product Description

1. LED Temperature display, resolution 0.1°C
2. Temperature control of the external system
3. PID temperature control
4. Low noise during the operation
5. Convenient to use on a table
6. Cooling coil (Option)
7. Adjustable high-temperature cut-off by heater and bath temperature sensor
8. Dry running protection
9. Protection from excess electric current
10. Water level sensor (Option)

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now