



WS-200(110V/60Hz)



\$1,371.00

Baths & Circulators

Circulating Water Bath
WS Series

Water Baths

The tank is using seamless stainless steel and the outside body is powder painting Screen-touch panel, micro-computer PID digital display Overheat protector The water bath with circulator using electric heat tube to unify temperature inside Add exterior circulator as requiring The sensor is PT-100Ω Gable cover (AC-100) sold separately for your option. The heater of stainless steel External Circulating Pump is Available Heating circulator sold separately for your option Optional Water level protector Optional Timer

[Add To Cart](#)

Product Description

The tank is using seamless stainless steel and the outside body is powder painting Screen-touch panel, micro-computer PID digital display Overheat protector The water bath with circulator using electric heat tube to unify temperature inside Add exterior circulator as requiring The sensor is PT-100Ω Gable cover (AC-100) sold separately for your option. The heater of stainless steel External Circulating Pump is Available Heating circulator sold separately for your option Optional Water level protector Optional Timer

Similar Products



Buy Now



Buy Now



Buy Now



Buy Now