



## WS-200(220V/50Hz)



\$1,371.00

Baths & Circulators

Circulating Water Bath WS Series

Water Baths

The tank is using seamless stainless steel and the outside body is powder painting Screen-touch panel, micro-computer PID digital display Overheat protector The water bath with circulator using electric heat tube to unify temperature inside Add exterior circulator as requiring The sensor is PT-100Ω Gable cover (AC-100) sold separately for your option. The heater of stainless steel External Circulating Pump is Available Heating circulator sold separately for your option Optional Water level protector Optional Timer

[Add To Cart](#)


### Product Description

The tank is using seamless stainless steel and the outside body is powder painting Screen-touch panel, micro-computer PID digital display Overheat protector The water bath with circulator using electric heat tube to unify temperature inside Add exterior circulator as requiring The sensor is PT-100Ω Gable cover (AC-100) sold separately for your option. The heater of stainless steel External Circulating Pump is Available Heating circulator sold separately for your option Optional Water level protector Optional Timer

### Similar Products



Buy Now



Buy Now



Buy Now



Buy Now